

Geobacillus stearothermophilus and Bacillus atrophaeus Dual Species Biological Indicator Spore Strip Product Comparison

This document refers to standard dual species spore strips of *Geobacillus stearothermophilus*, ATCC® 7953 and *Bacillus atrophaeus*, ATCC® 9372 NAMSA Product Code: STN-065R in comparison to Crosstex’s equivalent replacement Product Code: DS-100.


Product Code Reference and Shelf Life

| Product Description | NAMSA | | Crosstex | |
|--|--------------|------------|--------------|------------|
| | Product Code | Shelf Life | Product Code | Shelf Life |
| <i>Geobacillus stearothermophilus</i> 10 ⁵ Spore Strip | STN-065R | 30 months | DS-100 | 24 months |
| <i>Bacillus atrophaeus</i> 10 ⁶ Spore Strip | | | | |

Product Lot Comparison: NAMSA STN-065R versus Crosstex DS-100

| Product Characteristic | NAMSA | | | Crosstex | | |
|---|---------------------|---------------------|---------------------|----------------------|---------------------|---------------------|
| | Lot Number | | | Lot Number | | |
| | X912 | X937 | X938 | D35 | D36 | D38 |
| Population per strip | 1.9x10 ⁵ | 1.2x10 ⁵ | 1.1x10 ⁵ | 2.0x10 ⁵ | 2.2x10 ⁵ | 1.9x10 ⁵ |
| | 2.4x10 ⁶ | 1.8x10 ⁶ | 1.3x10 ⁶ | 1.9x10 ⁶ | 1.9x10 ⁶ | 1.8x10 ⁶ |
| Ethylene Oxide (EO) D value 600 mg/L, 54°C, 60% RH In Minutes | 3.3 | 3.4 | 3.2 | 3.0 | 3.0 | 2.8 |
| EO Survival Time In Minutes | 14.5 | 14.5 | 13.2 | Survives: 15 minutes | | |
| EO Kill Time In Minutes | 34.2 | 34.8 | 32.3 | Killed: 120 minutes | | |
| Dry Heat Survives 121°C | Not Certified | | | 30 minutes | | |
| Dry Heat D value at 160°C In Minutes | 1.9 | 1.9 | 1.6 | 2.3 | 2.3 | 2.1 |
| Dry Heat Survival Time at 160°C In Minutes | 8.4 | 8.1 | 6.6 | Not Certified | | |
| Dry Heat Kill Time at 160°C In Minutes | 19.7 | 19.4 | 16.1 | Killed: 120 minutes | | |

| Product Characteristic | NAMSA | | | Crosstex | | |
|---|---|------------------|------|----------------------|-----|-----|
| | Lot Number | | | Lot Number | | |
| | X912 | X937 | X938 | D35 | D36 | D38 |
| Dry Heat z value | Determined based on <i>D</i> values at 150°C, 160°C and 180°C | | | Not Certified | | |
| | 42°C | 42°C | 43°C | | | |
| Steam <i>D</i> value at 121°C In Minutes | 2.4 | 2.3 | 2.2 | 1.9 | 1.9 | 2.1 |
| Steam Survival Time at 121°C In Minutes | 7.9 | 7.1 | 6.7 | Survives: 5 minutes | | |
| Steam Kill Time at 121°C In Minutes | 22.2 | 20.8 | 19.8 | Killed: 15 minutes | | |
| Steam <i>D</i> value at 132°C In Minutes | Not Certified | Not Certified | 0.2 | Not Certified | | |
| Steam Survival Time at 132°C In Minutes | | | 0.7 | Survives: 20 seconds | | |
| Steam Kill Time at 132°C In Minutes | | | 1.8 | Killed: 2 minutes | | |
| Steam z value | Determined based on <i>D</i> values at 119°C, 121°C and 130°C | | | Not Certified | | |
| | 10°C | 10°C | 11°C | | | |



Crosstex is a FDA registered facility (registration number 1319130). Based on a product comparison and review of best practices, the Biological Indicator Spore Strips manufactured by Crosstex are equivalent to those manufactured by NAMSA Products Division and offer the following advantages:

- An improved peel area which allows for the Spore Strip to be more easily accessed and removed from the glassine packaging.
- Natural glassine material which does not contain dye.

| | NAMSA | Crosstex |
|---|--|--|
| Glassine Material | 30# Blue | 25# White |
| <p>Worldwide sources for glassine material are limited; thus, the material utilized by NAMSA and Crosstex-SPSmedical are the same with the exception of a blue dye. It is this dye which causes the glassine material NAMSA has utilized to appear blue, whereas Crosstex-SPSmedical utilizes the natural form of the material without the addition of the dye. Crosstex-SPSmedical switched from the dyed material to the natural white version several years ago without any negative impact to the product or its performance.</p> | | |
| Peel Area | Edges of glassine material align with indication to "PEEL HERE". | Edges of glassine material are offset, creating an obvious peel area which requires no indication. |
| Lot Number | 6-7 Digits printed on glassine | 4 Digits printed on glassine |
| Expiration Date | DD-MMM-YYYY | YYYY-MM-DD |
| Strip Dimension | 6 mm x 30 mm | 7 mm x 38 mm |
| Quantity | Boxes of 100 | Boxes of 100 |
| Population | 1.0 to 6.0 x 10 ^x | 1.0 to 5.0 x 10 ^x |
| EO D value | 2.5 to 5.8 minutes | 2.5 to 5.8 minutes |
| Dry Heat D value at 160°C | 1.0 to 3.0 minutes | 1.0 to 3.0 minutes |
| Dry Heat z value | ≥20°C | ≥20°C |

The ordering information associated with the products, during and post transition, will remain the same.

| Ordering Information | |
|----------------------|--|
| Telephone | 1 (800) 860-1888 |
| Fax | 1 (419) 666-1715 |
| Email | productorders@crosstex.com |

We understand that our clients have different needs based on your internal quality system requirements or Standard Operating Procedures (SOP's). Because of this, if there is something specific that is required, such as completion of specific forms or documents, samples, etc., please do not hesitate to let us know.

Thank you for being a valued client and we look forward to continuing to serve you.

If you have any questions, please do not hesitate to contact me.

Sincerely,



Julie Wheeler
Director Industrial Markets
jwheeler@crosstex.com
419-304-6535