

Culture Media - Manufacturing Transfer Summary

With the acquisition of NAMSA's Sterilization Monitoring Products Division, all manufacturing processes will be transferred to Crosstex. Culture Media has been manufactured by Crosstex for several years; thus, the validated processes and equipment currently being used to manufacture these products will continue to be utilized once the transition is complete.

Crosstex is a FDA registered facility (registration number 1319130). Based on a product comparison and review of best practices, the Culture Media manufactured by Crosstex is equivalent to that manufactured by NAMSA Products Division.

Product Comparisons

Culture Media for Use With *Bacillus atrophaeus* Biological Indicators and Inoculated Carriers

NAMSA Product Code: GMBTB-100

	NAMSA	Crosstex	
Product Codes	GMBTB-100	GMBTB-100	CM-100
Vial Dimensions	16 mm x 100 mm	13 mm x 100 mm	13 mm x 100 mm
Incubation Time	48 hours	48 hours	7 days
Media Transitioning (Positive)	Green to Yellow	Green to Yellow	Purple to Yellow
Packaging	100/box	100/box	100/box
Shelf Life	24 months	24 months	12 months

Culture Media for Use With *Geobacillus stearothermophilus* Biological Indicators and Inoculated Carriers

NAMSA Product Code: GMBCP-100

- NAMSA's Culture Media for use with *Geobacillus stearothermophilus* spores will be available under a new Crosstex product code.

	NAMSA	Crosstex
Product Codes	GMBCP-100	CM-100
Vial Dimensions	16 mm x 100 mm	13 mm x 100 mm
Incubation Time	24 hours	7 days
Media Transitioning (Positive)	Purple to Yellow	Purple to Yellow
Packaging	100/box	100/box
Shelf Life	24 months	12 months



PROTECTS®

(585) 359-0130 • Fax: (585) 359-0167
6789 West Henrietta Road • Rush, NY 14543 USA

Culture Media - Manufacturing Transfer Summary

The ordering information associated with the products, during and post transition will remain the same and as outlined below.

Ordering Information	
Telephone	1 (800) 860-1888
Fax	1 (419) 666-1715
Email	productorders@crosstex.com

As manufacturing activities begin to transfer to Crosstex, product may be branded under NAMSA or Crosstex and will be shipped from one of two locations depending on the branding of the product. NAMSA branded product will ship from the existing NAMSA facility (6750 Wales Road, Northwood, OH 43619) and Crosstex branded product will be shipped from the Rush, NY, facility (6789 W Henrietta Road, Rush, NY 14543).

Thank you for being a valued client and we look forward to continuing to serve you.